
Joypal

TOUCH PRO

HEATING MASSAGER

by SIDORA®

TOUCH PRO
HEATING MASSAGER

This manual should be well kept for easy reference.
It enables you to perform a safe and efficient operation.

Thank you for purchasing products. Before using the roduct, please read this manual thoroughly, particularly the Safety Precautions section, for the correct use of the product. This manual should be well kept for easy reference. It enables you to perform a safe and efficient operation.

CONTENTS

GUIDELINES FOR USE 4

BASIC PRODUCT INFORMATION 5

CHARGING / POWER INSTRUCTIONS 6

USE AND MAINTENANCE 6


HOW TO TROUBLESHOOT 7

WARNING 7


WARRANTY NOTES 9

CONTACT US 9


GUIDELINES FOR USE

PACKAGING INCLUDED:

Main Unit


Adapter

Instructions

GUIDELINES FOR USE

2000 mAh

Heating Button:

Press once to activate heating function. Press again to switch to low heat mode. Press third time to switch to high heat mode

Long press (3s) to turn on/off the power

Mode button:
Short press to switch modes

Mode 1: Comprehensive Alternating
Mode 2: Soothing Kneading Inward
Mode 3: Soothing Kneading Outward
Mode 4: Heating Mode

BASIC PRODUCT INFORMATION

Product Name:	TouchPro Heating Massager
Working Voltage:	11,1 V ⎓
Charging Voltage:	5 V ⎓
Operating Power:	15 W
Rated Input Power:	10 W
Battery Capacity:	2500 mAh
Charging Interface:	TYPE-C

Note: This product is a household product, please do not use it for commercial purposes, as we are constantly improving our products, so if there is any change ind the technical specifications, anger without notice. The picture of the manual is for reference only, the product is subject to the real thing.

CHARGING / POWER INSTRUCTIONS

Status of light	Description
Charging	The indicator light around the power button flashes slowly
Fully charged	The indicator light around the power button remains steadily on.
Low battery	The indicator lights on both buttons flash rapidly for 5 seconds and then turn off.

- This product can be used for 10 min in a single use, and it is recommended do use it no more than 2 times a day.
- This product can not be used while charging, please unplug the power adapter after charging.
- The product uses fast charging technology, when charging the hos and power adapter slightly warm, is a normal phenomenon.
- Please use our charging cable to charge the product and pay attention to the safety of electricity.
- Check the power supply, if the power supply voltage exceeds the specified range, the product cannot be charged.

USE AND MAINTENANCE

- Frequency of use, it is recommended to use the product no more than twice a day, and the time of each use should be controlled within 15 min.
- Before use, please do not operate the product when your hands are wet
- Product cleaning, if you need to clean the appearance, should be wiped with a soft dry cloth, do not rinse or wipe with water.
- Use/storage place: please keep the product away from water, fire source, strong sunlight, corrosive environment
- When the storage period exceeds 3 months, it should be opened and checked to recharge the machine before use.

HOW TO TROUBLESHOOT

Status	Reason	General Solutions
Unable to turn on	Whether the battery power is low	If yes, please change it in time
Massage strength is weak	Whether the battery power is weak	If yes, please change it in time
	Whether the device is in right place	Check whether the massage head is in a good fit with the top of the neck
Can not charging or no charging light	Whether the USB cable is connected to the main unit and the power adapter	If yes, please connect it in time
“Buzzing“ and “clattering” sound of the product	Sound made during normal operation of the product	No processing required

If the issue cannot be resolved after investigation, please contact online customer service on the purchasing platform. If there is a malfunction during operation, please perform the following checks. If the malfunction cannot be eliminated or the products till cannot operate normally, please cut off the power and have it repaired by professional technicians. Do not disassemble and assemble this product on your own.

WARNING

Please read the instructions carefully and do not use this product continuously for an extended period to prevent skin damage and muscle pain. Do not use the massager while driving Do not use if you have had incomplete recovery from surgery on the massage area or if you have skin damage, redness, swelling, inflammation, bruising or other discomfort. Do not use if you have heart disease, acute Illness, skin disease, malignant tumor, blood disease, osteoporosis, soft tissue disease, stroke, or other physical abnormalities.

If you have mobility, physical, or mental issues, consult your doctor before use. Pregnant women and individuals during menstruation are not recommended to use the massager without their doctor’s advice.

Children under 12 years of age are not recommended to use the massager. Children over 12 years of age must be supervised by an adult.

This product is designed for personal electronic health purposes and is not intended as a medical device. It should not for diagnostic, therapeutic, or commercial purposes. Before use this product, please remove any jewelry or other items that may obstruct the massage area to ensure safe and effective use. If you experience any pain or discomfort in the massage area during use, discontinue use immediately.In case of any abnormal behavior or malfunction of the product, please disconnect the power and contact our customer support team via the provided email address.

WARNING

Do not attempt to disassemble or repair the product yourself. Please be aware that this product has heat-generating areas on its surface.
Exercise caution when using it, especially if you are less sensitive to heat perception.

SPECIAL NOTE

We highly recommend that the following individuals exercise caution and seek guidance and advice from a healthcare professional before using this product.

1. Individuals suffering from malignant tumors

2. Patients with heart conditions or those using pacemakers

3. Individuals diagnosed with osteoporosis

4. Patients with recent physical trauma or skin diseases

5. Pregnant women, women in labor, or individuals experiencing their physiological periods

6. People with back problems, spinal abnormalities, or spinal twists

7. Individuals who are seriously physically or mentally handicapped or undergoing treatment for an illness

Safety Precautions for Specific Environments
When using this product in specific environments, please adhere to all relevant safety regulations to ensure your safety.

WARRANTY NOTED

RESPECTED USERS:

The warranty period for this product may vary according to the country of purchase, and we strive to maximize it within the boundaries set by local laws. As a standard, the warranty period for this product is 1 year from the date of purchase.
During the warranty period, if you encounter any quality issues with the product, we are committed to providing free repair or replacement services. Our company will bear the cost of replacing damaged parts or even replacing the entire machine, depending on the situation. Please note that damages resulting from improper use or irregular operation of the machine are not covered by the warranty. Ownership of the damaged product remains with our company unless a transfer of ownership is expressly agreed upon. Warranty or returns may not be available for the following circumstances:

1. .Product damage caused by unpredictable accidents such as fires, earthquakes, and floods.

2. Product failure or damage caused by negligence, misuse, or use under abnormal conditions.

3. Damage to the product caused by disassembly or improper storage, including damage caused by rodents or liquid infiltration.

4. Damage caused by the product being used in non-designated environments, including exceeding workload specifications, among other factors.

5. Product damage caused by circumstances beyond our control.

6. We are not liable for any incidental associated with the operation or use of the product, including loss of business profits, business interruptions, other damages resulting from product use, costs incurred for replacing equipment, devices, or services required for product use, or claims from your customers.

We are committed to adhering to applicable laws and regulations regarding liability for accidents, accidents in accordance with the relevant laws and regulations.

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Schlegel & Sidor International GmbH